





TGFB3 Antibody

Product Code	CSB-PA217147
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10600
Immunogen	Synthesized peptide derived from Human TGF β3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TGF β3 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	TRANSFORMING GROWTH FACTOR BETA-3;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFB3

Image



Western blot analysis of extracts from COS7 cells, using TGF β3 antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using TGF β 3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.