





TGFA Antibody

Product Code	CSB-PA799735
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01135
Immunogen	Synthesized peptide derived from Human TGF alpha.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGFA

Image



Western blot analysis of extracts from COS7 cells, using TGF alpha antibody.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma tissue using TGF alpha antibody.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.