





THRA Antibody

Product Code	CSB-PA593393
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10827
Immunogen	Synthesized peptide derived from Human thyroid hormone receptor α .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total thyroid hormone receptor α protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	C-erbA-alpha; EAR-7; EAR7; ERBA1; NR1A1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	THRA

Image



Western blot analysis of extracts from SKOV3 cells, using Thyroid Hormone Receptor α antibody.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma tissue using Thyroid Hormone Receptor α antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.