




TIMP3 Antibody

Product Code	CSB-PA218815
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35625
Immunogen	Synthesized peptide derived from internal of Human TIMP3.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TIMP3 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Metalloproteinase inhibitor 3 [Precursor]; TIMP-3; Tissue inhibitor of metalloproteinases 3; Protein MIG-5; TIMP3
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TIMP3

Image



Western blot analysis of extracts from HUVEC cells, using TIMP3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.