





NCOR1 Antibody

Product Code	CSB-PA035334
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75376
Immunogen	Synthesized peptide derived from Internal of Human NCoR1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total NCoR1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	KIAA1047; N-CoR1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NCOR1

Image



Western blot analysis of extracts from MDA-MB-435 cells, using NCoR1 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using NCoR1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.