





PIAS1 Antibody

Product Code	CSB-PA158229
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75925
Immunogen	Synthesized peptide derived from Internal of Human PIAS1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PIAS1 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AR interacting protein; DDXBP1; DEAD/H (Asp-Glu-Ala-Asp/His) box binding protein 1; DEAD/H box-binding protein 1; GBP
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIAS1

Image



Western blot analysis of extracts from MDA-MB-435 cells, using PIAS1 antibody.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using PIAS1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.