





PIAS2 Antibody

Product Code	CSB-PA546547
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75928
Immunogen	Synthesized peptide derived from Internal of Human PIAS2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total PIAS2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	ARIP3; Androgen receptor-interacting protein 3; DAB2-interacting protein; DIP; Miz1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PIAS2

Image



Western blot analysis of extracts from COS7 cells, using PIAS2 antibody.



Western blot analysis of extracts from HuvEc cells (Lane 2), JK cells (Lane 3) and cos-7 cells (Lane 4), using PIAS2 antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.