





# SENP8 Antibody

<b>Product Code</b>	CSB-PA117139
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96LD8
<b>Immunogen</b>	Synthesized peptide derived from Internal of Human SENP8.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Specificity</b>	The antibody detects endogenous levels of total SENP8 protein.
<b>Tested Applications</b>	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
<b>Form</b>	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Alias</b>	DEN1; NEDP; DENEDDYLASE 1; NEDD8-SPECIFIC PROTEASE 1; SENTRIN-SPECIFIC PROTEASE FAMILY
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SENP8

## Image



Western blot analysis of extracts from HuvEc cells, using SENP8 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue using SENP8 antibody.

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.