





TOP2B Antibody

Product Code	CSB-PA587958
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02880
Immunogen	Synthesized peptide derived from Internal of Human topoisomerase II β .
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total topoisomerase II β protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA topoisomerase II; beta isozyme; EC 5.99.1.3; TP2B; TOP2B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TOP2B


Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using Topoisomerase II β antibody.



Western blot analysis of extracts from Jurkat cells, using Topoisomerase II β antibody.



Western blot analysis of extracts from K562 cells (Lane 2), using TOP2B antibody. The lane on the left is treated with synthesized peptide.