





TOP2B Antibody

Product Code	CSB-PA587958
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02880
Immunogen	Synthesized peptide derived from Internal of Human topoisomerase II β .
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total topoisomerase II β protein.
Tested Applications	ELISA, WB, IHC; WB: 1:500-1:3000, IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA topoisomerase II; beta isozyme; EC 5.99.1.3; TP2B; TOP2B
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TOP2B


Image



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using Topoisomerase II β antibody.



Western blot analysis of extracts from Jurkat cells, using Topoisomerase II β antibody.



Western blot analysis of extracts from K562 cells (Lane 2), using TOP2B antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.