



Tyrosinase Antibody

Product Code	CSB-PA589649
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14679
Immunogen	Synthesized peptide derived from Internal of Human tyrosinase.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB; Recommended dilution: Western blotting: 1:500~1:3000
Form	Liquid
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TYR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.