





VEGFB Antibody

Product Code	CSB-PA959480
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49765
Immunogen	Synthesized peptide derived from Internal of Human VEGFB.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total VEGFB protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Vascular endothelial growth factor B; VEGF-B; VEGF-related factor; VRF; VASCULAR ENDOTHELIAL GROWTH FACTOR-RELATED FACTOR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VEGFB

Image



Western blot analysis of extracts from Jurkat cells, using VEGFB antibody.



Western blot analysis of extracts from K562 cells (Lane 2) and 293 cells (Lane 3), using VEGFB antibody. The lane on the left is treated with synthesized peptide.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.