




XRCC4 Antibody

Product Code	CSB-PA956465
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13426
Immunogen	Synthesized peptide derived from Internal of Human XRCC4.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total XRCC4 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA-repair protein XRCC4; X-ray repair cross-complementing protein 4;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	XRCC4

Image



Western blot analysis of extracts from COLO205 cells, using XRCC4 antibody.