




TNFRSF11B Antibody

Product Code	CSB-PA119257
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00300
Immunogen	Synthesized peptide derived from Internal of Human TR11B.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TR11B protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 11B [Precursor]; Osteoprotegerin; Osteoclastogenesis inhibitory factor; TNFRSF11B; OCIF
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF11B

Image



Western blot analysis of extracts from 293 cells, using TR11B antibody.