




TRIM38 Antibody

Product Code	CSB-PA275171
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00635
Immunogen	Synthesized peptide derived from internal of Human TRIM38.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TRIM38 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 38; RING finger protein 15; Zinc finger protein RoRet; TRIM38; RNF15
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TRIM38

Image



Western blot analysis of extracts from RAW264.7 cells and LOVO cells, using TRIM38 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.