





DLGAP1 Antibody

Product Code	CSB-PA778288
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14490
Immunogen	Synthesized peptide derived from internal of Human DLGP1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total DLGP1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DAP-1; DLGAP1; DLGP1; Disks large-associated protein 1; Guanylate kinase-associated protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DLGAP1

Image



Western blot analysis of extracts from HepG2 cells, using DLGP1 antibody.



Immunofluorescence analysis of HeLa cells, using DLGP1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.