




CDK2AP1 Antibody

Product Code	CSB-PA828207
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14519
Immunogen	Synthesized peptide derived from internal of Human CDKA1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total CDKA1protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Cyclin-dependent kinase 2-associated protein 1; CDK2-associated protein 1; Putative oral cancer suppressor; Deleted in oral cancer 1; DOC-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDK2AP1

Image



Western blot analysis of extracts from HeLa cells, using CDKA1 antibody.