




GAK Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA100194 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O14976 |
| Immunogen | Synthesized peptide derived from Internal of Human GAK. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse |
| Specificity | The antibody detects endogenous levels of total GAK protein. |
| Tested Applications | ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:200 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | Cyclin G-associated kinase; EC 2.7.11.1; kinase GAK; |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | GAK |

Image



Western blot analysis of extracts from 293 cells, using GAK antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.