





PDZD2 Antibody

Product Code	CSB-PA826920
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15018
Immunogen	Synthesized peptide derived from internal of Human PDZD2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PDZD2 protein.
Tested Applications	ELISA,IHC,IF;IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	PDZ domain-containing protein 2; PDZ domain-containing protein 3; Activated in prostate cancer protein; Processed PDZ domain-containing protein 2; PDZD2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PDZD2

Image



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue, using PDZD2 antibody.



Immunofluorescence analysis of HepG2 cells, using PDZD2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.