




EPHB6 Antibody

Product Code	CSB-PA017084
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15197
Immunogen	Synthesized peptide derived from internal of Human EPHB6.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total EPHB6 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EPB6; Ephrin type-B receptor 6 precursor; HEP; Tyrosine-protein kinase-defective receptor EPH-6;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPHB6


Image



Western blot analysis of extracts from HeLa cells, 293 cells and COLO cells, using EPHB6 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using EPHB6 antibody.



Immunofluorescence analysis of HeLa cells, using EPHB6 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.