




MAP3K13 Antibody

Product Code	CSB-PA029472
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43283
Immunogen	Synthesized peptide derived from internal of Human M3K13.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total M3K13 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 13; EC 2.7.11.25; Mixed lineage kinase; MLK; Leucine zipper-bearing kinase
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K13

Image



Western blot analysis of extracts from HeLa cells, HUVEC cells and HepG2 cells, using M3K13 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.