




ANP32C Antibody

Product Code	CSB-PA180795
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43423
Immunogen	Synthesized peptide derived from internal of Human AN32C.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total AN32C protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Acidic leucine-rich nuclear phosphoprotein 32 family member C; Phosphoprotein 32-related protein 1; Tumorigenic protein pp32r1; ANP32C; PP32R1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANP32C

Image



Western blot analysis of extracts from HUVEC cells, using AN32C antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.