





TNFSF12 Antibody

Product Code	CSB-PA152148
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43508
Immunogen	Synthesized peptide derived from internal of Human TNF12.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total TNF12 protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 12; TNF-related weak inducer of apoptosis; TWEAK; APO3 ligand; Tumor necrosis factor ligand superfamily member 12
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF12


Image



Western blot analysis of extracts from COLO cells, HUVEC cells and 293 cells, using TNF12 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using TNF12 antibody.



Immunofluorescence analysis of HeLa cells, using TNF12 antibody.