




TNFSF14 Antibody

Product Code	CSB-PA616434
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43557
Immunogen	Synthesized peptide derived from internal of Human TNF14.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TNF14 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 14; Herpesvirus entry mediator-ligand; HVEM-L; CD258 antigen; Tumor necrosis factor ligand superfamily member 14
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF14

Image



Immunofluorescence analysis of HeLa cells, using TNF14 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.