




AKAP10 Antibody

Product Code	CSB-PA250245
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43572
Immunogen	Synthesized peptide derived from Internal of Human AKAP10.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total AKAP10 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	A kinase anchor protein 10; mitochondrial [Precursor]; Protein kinase A-anchoring protein 10; PRKA10; Dual specificity A kinase-anchoring protein 2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP10

Image



Western blot analysis of extracts from mouse brain cells, using AKAP10 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.