





HR Antibody

Product Code	CSB-PA776465
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43593
Immunogen	Synthesized peptide derived from internal of Human HAIR.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total HAIR protein.
Tested Applications	ELISA,WB,IHC,IF;WB:1:500-1:3000,IHC:1:50-1:100,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Protein hairless; HR; HAIR;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HR


Image



Western blot analysis of extracts from 293 cells, using HAIR antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using HAIR antibody.



Immunofluorescence analysis of A549 cells, using HAIR antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.