





AIRE Antibody

Product Code	CSB-PA274868
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43918
Immunogen	Synthesized peptide derived from internal of Human AIRE.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total AIRE protein.
Tested Applications	ELISA, WB, IF; WB: 1:500-1:3000, IF: 1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Autoimmune regulator; Autoimmune polyendocrinopathy candidiasis ectodermal dystrophy protein; APECED protein; AIRE; APECED
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AIRE

Image



Western blot analysis of extracts from 293 cells, using AIRE antibody.



Immunofluorescence analysis of HeLa cells, using AIRE antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.