





ZEB2 Antibody

Product Code	CSB-PA031230
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60315
Immunogen	Synthesized peptide derived from internal of Human ZEB2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total ZEB2 protein.
Tested Applications	ELISA,WB,IHC;WB:1:500-1:3000,IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	SIP-1; SIP1; Smad interacting-protein 1; Smad-interacting 1; SMADIP1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZEB2

Image



Western blot analysis of extracts from HepG2 cells, using ZEB2 antibody.



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ZEB2 antibody.