




HBP1 Antibody

Product Code	CSB-PA097019
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60381
Immunogen	Synthesized peptide derived from Internal of Human HBP1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HBP1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HMG box-containing protein 1; HMG box-containing protein 1; High mobility group box transcription factor 1; HBP1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HBP1

Image



Western blot analysis of extracts from COS-7 cells, using HBP1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.