





SYNCRIP Antibody

Product Code	CSB-PA947601
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60506
Immunogen	Synthesized peptide derived from Internal of Human hnRNP Q.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total hnRNP Q protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoprotein Q; hnRNP Q; hnRNP-Q; Synaptotagmin-binding; cytoplasmic RNA-interacting protein
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SYNCRIP

Image



Western blot analysis of extracts from HUVEC cells, Jurkat cells and 293 cells, using hnRNP Q antibody.



Immunofluorescence analysis of HepG2 cells, using hnRNP Q antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.