




REV3L Antibody

Product Code	CSB-PA076411
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60673
Immunogen	Synthesized peptide derived from internal of Human DNA Polymerase ζ.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total DNA Polymerase ζ protein.
Tested Applications	ELISA,IHC;IHC:1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA polymerase zeta catalytic subunit; DPOLZ; DPOZ; EC 2.7.7.7; POLZ
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	REV3L

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using DNA Polymerase ζ antibody.