




REV3L Antibody

Product Code	CSB-PA076411
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60673
Immunogen	Synthesized peptide derived from internal of Human DNA Polymerase ζ .
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total DNA Polymerase ζ protein.
Tested Applications	ELISA, IHC; IHC: 1:50-1:100
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	DNA polymerase zeta catalytic subunit; DPOLZ; DPOZ; EC 2.7.7.7; POLZ
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	REV3L

Image



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using DNA Polymerase ζ antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.