




PFKFB2 Antibody

Product Code	CSB-PA796746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60825
Immunogen	Synthesized peptide derived from N-terminal of Human PFKFB2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PFKFB2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	6-phosphofructo-2-kinase; 6-phosphofructo-2-kinase/fructose-2; 6-biphosphatase 2; 6PF-2-K/Fru-2; 6-P2ase heart-type isozyme
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PFKFB2

Image



Western blot analysis of extracts from Jurkat cells, using PFKFB2 antibody.