




PFKFB2 Antibody

Product Code	CSB-PA796746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60825
Immunogen	Synthesized peptide derived from N-terminal of Human PFKFB2.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PFKFB2 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	6-phosphofructo-2-kinase; 6-phosphofructo-2-kinase/fructose-2; 6-biphosphatase 2; 6PF-2-K/Fru-2; 6-P2ase heart-type isozyme
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PFKFB2

Image



Western blot analysis of extracts from Jurkat cells, using PFKFB2 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.