





DKC1 Antibody

Product Code	CSB-PA820074
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60832
Immunogen	Synthesized peptide derived from Internal of Human Dyskerin.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Dyskerin protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	CBF5 homolog; DKC1; EC 5.4.99.-; Nucleolar protein NAP57;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	DKC1

Image



Western blot analysis of extracts from HeLa cells, using Dyskerin antibody.



Immunofluorescence analysis of HeLa cells, using Dyskerin antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.