




PFDN1 Antibody

Product Code	CSB-PA059337
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60925
Immunogen	Synthesized peptide derived from internal of Human PFDN1.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total PFDN1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Prefoldin subunit 1; PFDN1; PFD1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PFDN1

Image



Western blot analysis of extracts from RAW264.7 cells, using PFDN1 antibody.