




MAF Antibody

Product Code	CSB-PA796747
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75444
Immunogen	Synthesized peptide derived from Internal of Human Maf.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total Maf protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Transcription factor Maf; Proto-oncogene c-maf; MAF;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAF

Image



Western blot analysis of extracts from HuvEc cells, HepG2 cells, HeLa cells and COLO205 cells, using Maf antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.