



FOXH1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA791166 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75593 |
| Immunogen | Synthesized peptide derived from internal of Human FOXH1. |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Specificity | The antibody detects endogenous levels of total FOXH1 protein. |
| Tested Applications | ELISA, IF; IF: 1:100-1:500 |
| Form | Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Alias | Forkhead box protein H1; Forkhead activin signal transducer 1; Fast-1; hFAST-1; Forkhead activin signal transducer 2 |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | FOXH1 |

Image



Immunofluorescence analysis of A549 cells, using FOXH1 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.