



AKAP3 Antibody

Product Code	CSB-PA588403
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75969
Immunogen	Synthesized peptide derived from internal of Human AKAP3.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total AKAP3 protein.
Tested Applications	ELISA,IF;IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	AKAP110; FSP95; PRKA3; SOB1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP3

Image



Immunofluorescence analysis of NIH/3T3 cells, using AKAP3 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.