




FARP2 Antibody

Product Code	CSB-PA977826
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94887
Immunogen	Synthesized peptide derived from Internal of Human FIR.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total FIR protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	FERM; RhoGEF and pleckstrin domain protein 2; FGD1-related Cdc42-GEF; FIR; FRG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FARP2

Image



Western blot analysis of extracts from HuvEc cells, using FIR antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.