





HEXIM1 Antibody

Product Code	CSB-PA556810
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94992
Immunogen	Synthesized peptide derived from Internal of Human HES7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HES7 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Protein HEXIM1; Hexamethylene bis-acetamide-inducible protein 1; Cardiac lineage protein 1; Estrogen down-regulated gene 1 protein; Menage a quatre protein 1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HEXIM1

Image



Western blot analysis of extracts from NIH-3T3 cells, using HES7 antibody.



Immunofluorescence analysis of HeLa cells, using HES7 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.