





TNFSF15 Antibody

Product Code	CSB-PA983205
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95150
Immunogen	Synthesized peptide derived from C-terminal of Human TNFSF15.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total TNFSF15 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 15; Vascular endothelial cell growth inhibitor; TNF ligand-related molecule 1; TNFSF15; TL1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF15

Image



Western blot analysis of extracts from COLO205 cells, using TNFSF15 antibody.



Immunofluorescence analysis of HeLa cells, using TNFSF15 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.