





ZRANB2 Antibody

Product Code	CSB-PA104679
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95218
Immunogen	Synthesized peptide derived from internal of Human ZNF265.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total ZNF265 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Z265; ZIS; ZN265; Zinc finger protein 265; Zinc finger
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ZRANB2

Image



Western blot analysis of extracts from HeLa cells, using ZNF265 antibody.



Immunofluorescence analysis of HUVEC cells, using ZNF265 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.