





KAT7 Antibody

Product Code	CSB-PA980473
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95251
Immunogen	Synthesized peptide derived from internal of Human MYST2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total MYST2 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.3.1.48; HBO1; HBOA; MYST histone acetyltransferase 2;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KAT7

Image



Western blot analysis of extracts from Jurkat cells, using MYST2 antibody.



Immunofluorescence analysis of HUVEC cells, using MYST2 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.