




YEATS4 Antibody

Product Code	CSB-PA105916
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95619
Immunogen	Synthesized peptide derived from internal of Human GAS41.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Specificity	The antibody detects endogenous levels of total GAS41 protein.
Tested Applications	ELISA, WB; WB: 1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	YEATS domain-containing protein 4; Glioma-amplified sequence 41; Gas41; NuMA-binding protein 1; NuBI-1
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	YEATS4

Image



Western blot analysis of extracts from mouse lung cells, using GAS41 antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.