





MAP4K4 Antibody

Product Code	CSB-PA781476
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95819
Immunogen	Synthesized peptide derived from internal of Human MEKKK 4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous levels of total MEKKK 4 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EC 2.7.11.1; HPK/GCK-like kinase HGK; M4K4; MAP4K4; MAPK/ERK kinase kinase kinase 4
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MAP4K4

Image



Western blot analysis of extracts from COLO cells and HUVEC cells, using MEKKK 4 antibody.



Immunofluorescence analysis of A549 cells, using MEKKK 4 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.