





PTTG1 Antibody

Product Code	CSB-PA443065
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95997
Immunogen	Synthesized peptide derived from internal of Human PTTG1.
Raised In	Rabbit
Species Reactivity	Human
Specificity	The antibody detects endogenous levels of total PTTG1 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EAP1; PTTG; PTTG1; TUTR1;
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PTTG1

Image



Western blot analysis of extracts from 293 cells, Jurkat cells and HepG2 cells, using PTTG1 antibody.



Immunofluorescence analysis of HUVEC cells, using PTTG1 antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.