





NGFR Antibody

Product Code	CSB-PA038076
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08138
Immunogen	Synthesized peptide derived from internal of Human TNR16.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TNR16 protein.
Tested Applications	ELISA,WB,IF;WB:1:500-1:3000,IF:1:100-1:500
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	Gp80-LNGFR; NGF receptor; TNFRSF16; TNR16; low affinity neurotrophin receptor p75NTR
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	NGFR

Image



Western blot analysis of extracts from HepG2 cells, using TNR16 antibody.



Immunofluorescence analysis of A549 cells,
using TNR16 antibody.