




HMGB1 Antibody

Product Code	CSB-PA049959
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09429
Immunogen	Synthesized peptide derived from internal of Human HMGB1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total HMGB1 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	HMG-1; HMG1; High mobility group protein 1; SBP-1; amphoterin
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HMGB1

Image



Western blot analysis of extracts from Jurkat cells, HepG2 cells, COLO205 cells and 293 cells, using HMGB1 antibody.