




TACSTD2 Antibody

Product Code	CSB-PA170174
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09758
Immunogen	Synthesized peptide derived from internal of Human TACD2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous levels of total TACD2 protein.
Tested Applications	ELISA,WB;WB:1:500-1:3000
Form	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Alias	EGP-1; GA733-1; M1S1; TACD2; TROP2
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TACSTD2

Image



Western blot analysis of extracts from A549 cells and NIH-3T3 cells, using TACD2 antibody.